



LIFEVAC
SOUTH AFRICA



WHEN FIRST AID AND ALL OTHER
METHODS FAIL IN A CHOKING EMERGENCY

LIFEVAC SAVES LIVES

WHY CHOOSE LIFEVAC

TO PROTECT THE PEOPLE YOU CARE FOR IN A CHOKING EMERGENCY?

LIFEVAC IS **SIMPLE, EFFECTIVE, NON-INVASIVE** AND MOST IMPORTANTLY **SAFE TO USE.**

LifeVac are the leaders in airway clearance device technology, being the only ACD with peer reviewed medical publications, independent medical testing and medical abstracts proving safety, effectiveness, usability and lives saved.

WHAT IS LIFEVAC?

LifeVac is a self-powered, single patient portable suction device developed for clearing the upper airway, intended to be used in an emergency when standard current choking protocol has been followed without success.

- 1 Interchangeable sized masks to fit a casualties facial features, as one size does not fit all.
- 2 One-way valve prevents any air being expelled through interchangeable sized masks.
- 3 Translucent bellows, makes it easy to identify if the obstruction enters this area.
- 4 Easy to hold handle for secure grip.


THE NEGATIVE PRESSURE GENERATED BY THE FORCE OF THE SUCTION IS 3 TIMES GREATER THAN THE HIGHEST RECORDED FIRST-AID MEASURE.

LifeVac generates over 300 mmHg of suction. The duration of suction is minimal so LifeVac is safe and effective. No prescription is required. LifeVac requires no maintenance.




EASY AS


1 PLACE



2 PUSH



3 PULL





ONLY AIRWAY CLEARANCE DEVICE WITH INDEPENDENT MEDICAL TESTING

Peer reviewed medical publications, peer reviewed abstracts proving safety, effectiveness and lives saved.

KEY FEATURES

- LifeVac is a non-invasive, portable airway clearance device.
- Interchangeable sized masks, clearly identified by colour coded rings.
- No risk of pushing the tongue or obstruction back in a panic situation.
- No risk of oral damage.
- Generates over 326mm Hg of suction, safely and effectively dislodging the obstruction.
- Can be used for full and partial obstructions.
- Saved many lives around the world from choking to death.
- Can be applied if someone is standing, sitting or laying down.

AVAILABLE OPTIONS

- LifeVac comes in three different variations, Standard Kit, Wall mounted Kit and LifeVac EMS Kit
- Different Mask Sizes are available to purchase separately as well.

CONTACT

011 450 0096

WWW.LIFEVAC.CO.ZA

LIFEVAC@VENTRIA.CO.ZA



STANDARD LIFEVAC KIT INCLUDES:

- Toddler
- Small Adult/Child
- Large Adult Masks



WALL MOUNTED KIT INCLUDES:

- Toddler
- Small Adult/Child
- Medium Adult
- Large Adult Masks

EMS KIT INCLUDES:

- Toddler
- Small Adult/Child
- Medium Adult
- Large Adult
- Extra Large Adult Masks



LIFEVAC IS FDA, MHRA REGISTERED AND CE ACCREDITED